

# Certificate of Analysis

(Attached)

## **SOURCE MATERIAL TESTED AND FOUND NEGATIVE FOR ANTIBODY TO HIV AND ANTIBODY FOR HBV AND FOR HCV AND FOR HTLV**

### **Caution**

1. These pathogen negative materials and derivative products are potential infectious materials.
2. There is a possibility that testing result for HBV, HCV, HIV, and HTLV may be false negative.
3. Current clinical testing kit may not detect mutated version of these virus.
4. Products may contain some other unknown pathogens.

### **Please handle this material or product as potential infectious**

1. Inform any personnel, customers, and users, who contact these materials or products, that these pathogen negative materials and derivative products are potential infectious materials.
2. Adopt sterilizing procedure to handle pathogen negative materials and derivative products whenever and wherever possible.
3. Protect the personnel who contact these pathogen negative materials and derivative products, and protect environment during using these materials.
4. Treat pathogen negative materials and derivative products as potential infectious materials during experiment, testing, manufacture, production, packaging, shipping, handling, using, and disposing.
5. Treat gel buffer, tips, tubes, bench top, containers, etc. contaminated by these materials or products with sufficient 2% bleach for 15 to 60 minutes to eliminate all potential infection