

Data Sheet

Product Name: cDNA Library

Catalog No.: B1234035

Lot No.: A700101

Species: Human Mouse Rat Monkey Plant Cell line Dog

Tissue Type: Normal Adult Fetal Tumor Disease

Tissue Name: Brain

Vector: BioExpress

Cells used for Electroporation: DH10B *E. coli* bacterial cells, T1 phage resistant

Diluent: 80% SOC:20% Glycerol

Titer: 5 x10⁹ cfu/ml

Average insert size: 2.1 kb

Percent recombinants: 100%

Primer(s) for sequencing: M13FWD 5'-GTA AAA CGA CGG CCA GT-3'
M13REV 5'-CAG GAA ACA GCT ATG AC-3'

Donor Information:

Male: 23-27 year(s) old

Pool of 3 donor(s)

Female: _____ year(s) old

Pool of _____ donor(s)

Pathological Diagnosis:

Components

1. 1 ml library
2. cDNA library user manual
3. Certificate of Analysis
4. Data Sheet

FOR RESEARCH USE ONLY

APPROVED BY: _____